


A-A (0.09 : 1)



NOTES :

CUSTOMER :

- SUMIT DRAWING FOR APPROVAL
- 1.0 CONSTRUCTION: EPOXY COATES MILD STEEL CABINET, CART, POLYCARB SIDES PANELS AND HEAD UNIT.
- 2.0 FRONT SLIDING SASH FUME HOOD.
- 3.0 CENTRIFUGAL FAN.
- 4.0 POWER SUPPLY 110 V 60 HZ.
- 5.0 MONITORING FLOW AIR FLOW.
- 6.0 ELECTRONIC ALARM SYSTEM .
- 7.0 HEAVY DUTY CART WITH TWO FRONT LOCKABLE WHEELS.
- 8.0 INCLUDED:
- 8.1 ELECTRONIC FILTER SATURATION ALARM.
- 8.2 PURAIR ACI MAIN FILTER.
- 8.3 PURAIR 20 SPILLAGE TRAY.
- 9.0 W/O

Designed by CFG	Checked by	Approved by	Date 1/14/2014	Date
			ECO-20 UNIT	
			ECO20.GAssem	Sheet 1 / 2