



ELECTRICAL CONNECTION

MAIN FILTERS

DUPLEX ELECTRICAL OUTLET

BACKUP FILTERS

GAS VALVE

FUME HOOD CONTROLLER NETWORK CABLE

WATER NECK FAUCET

NOTES :

- 1.0 CUSTOMER
- 2.0 SUBMITTED DRAWING FOR APPROVAL.
- 3.0 STANDARD FEATURES INCLUDED.
- 3.1 CONSTRUCTION, EPOXY COATED MILD STEEL.
- 3.2 FRONT SASH AND SIDES PANELS MADE IN CLEAR TEMPERED GLASS.
- 3.3 POWER SUPPLY 110 V 60 HZ.
- 3.5 ON / OFF SWICHT ON CONTROL PANEL.
- 3.6 STANDARD ACCESORIES.

323,00 / 12.72"

965,20 / 38.0"

1244,60 / 49.0"

1297,94 / 51.10" Aprox.

1498,60 / 59.0"

CUP SINK

GAS VALVE FITTING

SLIDING FRONT SASH

Designed by Carlos Granados	Checked by	Approved by	Date	Date 3/11/2010
			ECO25 + SINK	
			ECO25-Gasemt	Sheet 1 / 1