



- NOTES :**
- 1.0 CUSTOMER
 - 2.0 SUBMIT DRAWING FOR APPROVAL.
 - 3.0 STANDARD FEATURES INCLUDED.
 - 3.1 CONSTRUCTION, EPOXY COATED MILD STEEL.
 - 3.2 FRONT SASH AND SIDES PANELS MADE IN CLEAR TEMPERED GLASS.
 - 3.3 POWER SUPPLY 110 V 60 HZ.
 - 3.5 OPTIONAL MONITAIR.
 - 4.0 PAINT : YUKX WHITE PPL 51273
 BLACK ONYX COLOR EPOXY RESIN COUNTER UP .

Designed by CFG	Checked by	Approved by	Date 6/1/2011	Date
ECO30 UNIT			Edition	Sheet 1 / 1
			ECO30.GAssem	