



- NOTES:**
- STANDARD ON ALL UNITS:
 - 1.0 CONSTRUCTION: EPOXY COATED MILD STEEL HEAD UNIT AND COLUMNS.
 - 2.0 FRONT, SIDES AND BACK PANELS MADE ON CLEAR ACRYLIC OF 1/4".
 - 3.0 CENTRIFUGAL FAN
 - 4.0 POWER SUPPLY 110V 60Hz
 - 5.0 SWITCHES, MAIN ON/OFF
 - 6.0 MONITORING, LOW AIRFLOW

| | | | | |
|--------------------|------------|-------------|----------------------|----------------|
| Designed by CFG | Checked by | Approved by | Date 1/30/2013 | Date |
| | | | P5.24SL.XT UNIT | |
| | | | P5.24SLXCTR12GIAssem | Sheet 1 / 1 |