



- NOTES:**
- A.0 STANDARD ON ALL UNITS:
  - 1.0 CONSTRUCTION: EPOXY COATED MILD STEEL HEAD UNIT AND COLUMNS.
  - 2.0 FRONT, SIDES AND BACK PANELS MADE OF 1/4" CLEAR ACRYLIC.
  - 3.0 COLOR EPOXY COATED BLUE
  - 4.0 CENTRIFUGAL FAN.
  - 5.0 POWER SUPPLY:
  - 5.1 FOR 120 VOLT 60 HZ.
  - 5.2 FOR 230 VOLT 50 HZ.
  - 6.0 SWITCHES, MAIN ON/OFF.
  - 7.0 MONITORING, LOW AIRFLOW.
  - 8.0 HINGED FRONT SASH.
  - 9.0 ELECTROSTATIC PREFILTER.
  - 10.0 TWO PASS THROUGH PORTS.
  - 11.0 FILTER DOOR WITH KEY.
  - 12.0 MAIN FILTER IS APPLICATION-SPECIFIC; MUST ORDER SEPARATELY.

Designed by CFG	Checked by	Approved by	Date 3/20/2015	Date
P536XT-R15 METAL BAFFLE			Edition	Sheet 1 / 1
			P536XTR15-MB	