



**NOTES:**

- A.0 STANDARD ON ALL UNITS:
- 1.0 CONSTRUCTION: EPOXY COATED MILD STEEL HEAD UNIT AND COLUMNS.
- 2.0 FRONT, SIDES AND BACK PANELS MADE OF 1/4" CLEAR ACRYLIC.
- 3.0 COLOR EPOXY COATED BLUE
- 4.0 CENTRIFUGAL FAN.
- 5.0 POWER SUPPLY:
- 5.1 FOR 120 VOLT 60 HZ.
- 5.2 FOR 230 VOLT 50 HZ.
- 6.0 SWITCHES, MAIN ON/OFF.
- 7.0 MONITORING, LOW AIRFLOW.
- 8.0 HINGED FRONT SASH.
- 9.0 ELECTROSTATIC PREFILTER.
- 10.0 TWO PASS THROUGH PORTS.
- 11.0 FILTER DOOR WITH KEY.
- 12.0 MAIN FILTER IS APPLICATION-SPECIFIC; MUST ORDER SEPARATELY.

Designed by CFG	Checked by	Approved by	Date 3/20/2015	Date
P536XT-R15 METAL BAFFLE				
P536XTR15-MB			Edition	Sheet 1 / 1