




8		7		6	
REVISION HISTORY					
ZONE	REV	DESCRIPTION	DATE	APPROVED	
WORK SURFACE	1	SS WORK SURFACE 40 MM HIGH	1/20/2013	P CH	
TOP FRONT	2	CONTROL OANEL POSITION	7/25/2018	P CH	



A-A (1 : 16)



NOTES :

- 1.0 CUSTOMER
- 2.0 SUBMIT DRAWING FOR APPROVAL.
- 3.0 STANDARD FEATURES INCLUDED.
 - 3.1 CONSTRUCTION, EPOXY MICROgone
ANTIMICROBIAL COATING MILD STEEL.
 - 3.2 STAINLESS STEEL WORK SURFACE .
 - 3.3 SIDES PANELS MADE IN TEMPERED GLASSES.
 - 3.4 FANFIL-FARR ULPA MAIN FILTER.
 - 3.5 INTERNAL LIGHTNING.
 - 3.6 INTERNAL GFCI OUTLET.
 - 3.7 FAN SPEED CONTROL.
 - 3.8 OPTIONAL; BASE STAND OR MOVIL CART.
 - 3.9 CONTROL PANEL: FAN'S ON/OFF, MINIHELIC
ULPA PRESSURE GAUGE, GFCI OUTLET, AND
LIGHTS SWITCHES.
- 4.0 EBMFAPST ROTOR BLOWER.
- 5.0 POWER INLET 120 V 60 HZ.
- 6.0 OPTIONAL; HANGING BAR WITH "S" HOOKS.
- 7.0 W/O

Designed by CFG	Checked by	Approved by	Date	Date 4/18/2011	
			VLF-72 VERTICAL LAMINAR FLOW		
			VLF-72.GAssem	Edition 2	Sheet 1 / 1