




- NOTES:**
- A.0 STANDARD ON ALL UNITS:
 - 1.0 CONSTRUCTION: EPOXY COATED MILD STEEL HEAD UNIT AND FRAME.
 - 2.0 FRONT, SIDES AND BACK PANELS MADE OF 1/4" CLEAR POLYCARBONATE.
 - 3.0 MAIN SWITCH ON/OFF.
 - 4.0 CENTRIFUGAL FAN.
 - 5.0 POWER SUPPLY:
 - 5.1 FOR 120 VOLT 60 HZ.
 - 5.2 FOR 230 VOLT 50 HZ.
 - B.0 MAIN FILTER IS APPLICATION-SPECIFIC; MUST ORDER SEPARATELY.

Designed by	Checked by	Approved by	Date	Date	
			6/25/2019		
			MINI SINGLE SAFESTORE CSC		
			Mini-Stg01.GAssy	Edition	Sheet
					1 / 1