



FILTER ACCESS DOOR WITH VIEWER

POWER INLET

A-A (1 : 6)

ELECTRIC PORTS (2)

MAIN ON/OFF SWITCH

POLYPRO SPILLAGE TRAY (OPTIONAL)

NOTES:

- STANDARD ON ALL UNITS:
- 1.0 CONSTRUCTION: EPOXY COATED MILD STEEL HEAD UNIT AND COLUMNS.
- 2.0 FRONT, SIDES AND BACK PANELS MADE ON CLEAR ACRYLIC OF 1/4".
- 3.0 CENTRIFUGAL FAN
- 4.0 POWER SUPPLY 110V 60Hz
- 5.0 SWITCHES, MAIN ON/OFF
- 6.0 MONITORING, LOW AIRFLOW

Designed by CFG	Checked by	Approved by	Date	Date 1/29/2013
			P5.24XT UNIT	
			P5.24AXTR12-GAssem	Sheet 1 / 1