



- NOTES :
- STANDARD ON ALL UNITS
- 1.0 CONTRUCTION: EPOXY COATED MILD STEEL CABINET, CART, SIDE PANELS AND HEAD UNIT.
  - 2.0 FRONT, SIDES, BACK, DIFUSER AND TOP PANELS MADE IN CLEAR ACRYLIC OF 1/4".
  - 3.0 CENTRIFUGAL FAN.
  - 4.0 POWER SUPPLY 110 V 60 HZ.
  - 5.0 RATING ON MOTOR 200 CFM .
  - 6.0 FACE VELOCITY 60-100 LFM.
  - 7.0 PRE-FILTER 99.6% EFFICIENCY.
  - 8.0 EDU CARBON FILTER.
  - 9.0 MONITORING LOW AIR FLOW.
  - 10.0 WATER & GAS SERVICES, OPTIONALS.
  - 11.0 HEAVY DUTY CART WITH TWO LOCABLE WHEELS.

Designed by CFG	Checked by	Approved by	Date 8/31/2009	
EDU WORK STATION UNIT				
EDU.Mobile.Assem(Monitair)			Edition R 01	Sheet 1 / 1