


REVISION HISTORY				
ZONE	REV	DESCRIPTION	DATE	APPROVED
CART	1	CUT OUT SINK	8/29/2017	P CH



- NOTES :
- STANDARD ON ALL UNITS
 - 1.0 CONSTRUCTION: EPOXY COATES MID STEEL CABINET, CART, SIDE PANELS AND HEAD UNIT.
 - 2.0 FRONT, SIDES, BACK, DIFUSER AND TOP PANELS MADE IN CLEAR ACRYLIC OF 1/4".
 - 3.0 CENTRIFUGAL FAN.
 - 4.0 POWER SUPPLY 110 V 60 HZ.
 - 5.0 RATING ON MOTOR 200 CFM .
 - 6.0 FACE VELOCITY 60-100 FPM.
 - 7.0 PRE-FILTER 99.6% EFFICIENCY.
 - 8.0 EDU CARBON FILTER.
 - 9.0 MONITORING LOW AIR FLOW.
 - 10.0 HEAVY DUTY CART WITH TWO LOCABLE WHEELS.

Designed by CFG	Checked by	Approved by	Date 12/28/2015	
			EDU M-60 WORK STATION UNIT	
EDU60.WorkS.Assem			Edition 1	Sheet 1 / 1