




REVISION HISTORY				
ZONE	REV	DESCRIPTION	DATE	APPROVED
SIDE PANELS	1	TEMPERED GLASSES	3/22/2015	P CH




B-B (1 : 10)



A-A (1 : 10)

NOTES:

- STANDARD ON ALL UNITS:
- 1.0 CONSTRUCTION: EPOXY COATED MID STEEL SIDES PANELS AND HEAD UNIT.
 - 2.0 FRONT AND BACK PANELS MADE ON CLEAR ACRYLIC OF 1/4".
 - 3.0 TEMPERED GLASSES ON SIDE PANELS.
 - 4.0 CENTRIFUGAL FAN
 - 5.0 POWER SUPPLY 110V 60Hz
 - 6.0 MAIN SWITCH ON/OFF
 - 7.0 MONITORING, LOW AIRFLOW



Designed by CFG	Checked by	Approved by	Date 12/10/2008	
			PURAIR 10XL DUCTLESS FUME HOOD	
			Purair10XL-UL.Assem	Edition 1
				Sheet 1 / 1