




NOTES:

STANDARD IN ALL UNITS:

- 1.0 CONSTRUCTION: EPOXY COATED MID STEEL SIDES PANELS AND HEAD UNIT.
- 2.0 FRONT, SIDES AND BACK PANELS MADE ON CLEAR ACRYLIC OF 1/4".
- 3.0 TEMPERED GLASSES ON SIDE PANELS.
- 4.0 CENTRIFUGAL FAN.
- 5.0 POWER SUPPLY 110V 60Hz.
- 6.0 MAIN SWITCH ON/OFF.
- 7.0 MONITORING, LOW AIRFLOW.

Designed by CFG	Checked by	Approved by	Date 3/6/2014	
			PURAIR 15 DUCTLESS FUME HOOD	
			Purair15-UL.Assem	Edition 1 Sheet 1 / 1