




REVISION HISTORY				
ZONE	REV	DESCRIPTION	DATE	APPROVED
MF LEVEL	1	LIGHTS BASE POSITION	3/29/2017	P CH




A-A ( 1 : 10 )



B-B ( 1 : 10 )

- NOTES:
- STANDARD ON ALL UNITS:
- 1.0 CONSTRUCTION:  
ENCLOSURE BUILT 3/8" POLYPRO SEAMLESS, NON-POROUS, WATER-TIGHT INSIDE,THERMALY FUSE. FAN MODULE BUILT EPOXY COATED MILD STEEL.
  - 2.0 FRONT AND SIDES PANELS MADE IN CLEAR POLYCARB 1/4".
  - 3.0 CENTRIFUGAL FAN
  - 4.0 POWER SUPPLY 110V 60Hz
  - 5.0 SWITCHES, MAIN ON/OFF
  - 6.0 MONITORING, LOW AIRFLOW



Designed by cfg	Checked by	Approved by	Date	Date 10/25/2017	
			POLYPRO PCR 36		
			PPCR36-Gassembly	Edition 1	Sheet 1 / 1