



| REVISION HISTORY |     |                            |           |          |
|------------------|-----|----------------------------|-----------|----------|
| ZONE             | REV | DESCRIPTION                | DATE      | APPROVED |
| WORK SURFACE     | 1   | SS WORK SURFACE 40 MM HIGH | 1/10/2013 | P CH     |
| FRONT TOP        | 2   | CONTROL PANEL POSITION     | 2/20/2018 | P CH     |

F  
E  
D  
C  
B  
A



- NOTES :
- 1.0 CUSTOMER
  - 2.0 SUBMIT DRAWING FOR APPROVAL.
  - 3.0 STANDARD FEATURES INCLUDED.
  - 3.1 CONSTRUCTION, EPOXY MICROgone ANTIMICROBIAL COATING MILD STEEL.
  - 3.2 STAINLESS STEEL WORK SURFACE .
  - 3.3 SIDES PANELS MADE IN TEMPERED GLASSES.
  - 3.4 CANFIL-FARR ULPA MAIN FILTER.
  - 3.5 INTERNAL LIGHTNING.
  - 3.6 INTERNAL GFCI OUTLET.
  - 3.7 FAN SPEED CONTROL.
  - 3.8 OPTIONAL; BASE STAND OR MOVIL CART.
  - 3.9 CONTROL PANEL: FAN'S ON/OFF, MINIHILIC ULPA PRESSURE GAUGE, GFCI OUTLET, AND LIGHTS SWITCHES.
  - 4.0 EBMPAPST ROTOR BLOWER.
  - 5.0 POWER INLET 120 V 60 HZ.
  - 6.0 OPTIONAL; HANGING BAR WITH "S" HOOKS.
  - 7.0 HEPA FILTER 915 X 610 X 66.



A-A ( 1 : 14 )

B-B ( 1 : 14 )

|                    |            |             |                              |                          |
|--------------------|------------|-------------|------------------------------|--------------------------|
| Designed by<br>CFG | Checked by | Approved by | Date<br>2/9/2012             |                          |
|                    |            |             | VLF-36 VERTICAL LAMINAR FLOW |                          |
|                    |            |             | VLF-36.GAssem                | Edition 2<br>Sheet 1 / 1 |

F  
E  
D  
C  
B  
A