


REVISION HISTORY				
ZONE	REV	DESCRIPTION	DATE	APPROVED
FAN LEVEL	1	UL FAN MODULE	7/06/2010	P CH
MAIN FILTER	2	METAL CLAMPS	8/12/2011	P CH
RIGHT SIDE	3	METAL SHOEBOX	3/18/2012	P CH
TOP COVER	4	UL FAN METAL MODULE	11/8/2017	P CH



THREADED BRONZE
1/2" Ø STD FEMALE
WATER HOSE
INLET CONECTOR

DRAINAGE 1 1/2" Ø
PVC PIPE
AND FITTINGS



MAIN CONTROL
WATER VALVE

THREADED BRONZE
1/2" Ø STD FEMALE
DRAIN HOSE
OUTLET CONECTOR


WASTED
PUMP

HOSE WATER
OUTLET TO
SPRAY VALVE

B-B (1 : 10)



A-A (1 : 10)



NOTES:

1.0 CUSTOMER:

2.0 SUBMIT DRAWINGS FOR APPROVAL.

3.0 INCLUDED:

3.1 STANDARD ITEMS APPLIED.

3.2 PLUMBING PACKAGE.

3.3 SPRAY VALVE AND 1/2" X 6' HOSES

FOR DRAINAGE AND COLD WATER

WITH BRONZE FEMALE CONECTORS.

4.0 CABINET BUILT WITH 3/8" WHITE

POLYPROPYLENE, SEAMLESS, NON-POROUS

WATERTIGHT INSIDE, THERMALLY FUSE.

5.0 STANDARD MOLDED PLUG NEMA 4-14P

6.0 ELECTRICAL REQUIREMENTS 110V 6A 60 HZ

Designed by CFG	Checked by	Approved by	Date 10/17/2011	
DC36XT.2000UL SAFEKEEPER				
DC36XTUL-GAssem			Edition 4	Sheet 1 / 1